

***DAR-9400f/DAR-9500f***  
**with Gateway Option**  
**DICOM Conformance Statement**

Revision E

Revision Date: April, 2015



| Date      | Rev.  | Comment   |
|-----------|-------|---|
| 2008.1.15 | First | First Revision  |
| 2008.6.25 | A     | Revise Chapter 6  |
| 2012.3.16 | B     | Apply to DAR-9500   |
| 2013.9.13 | C     | Apply to RDSR   |
| 2014.4.4  | D     | Apply to DAR-9400f Ver.4.2 and DAR-9500f Ver.5.2.40<br>Add description of usage to MPPS and RDSR.<br>Add information object definition of RDSR. |
| 2015.4.23 | E     | Apply to DAR-9500f Ver.6.0<br>Add private tag for StentView and 3D application.   |

## TABLE OF CONTENTS

|   |           |
|---|-----------|
| <b>1. INTRODUCTION</b> .....  | <b>1</b>  |
| 1.1. PURPOSE OF THIS DOCUMENT .....   | 1         |
| 1.2. SOURCES FOR THIS DOCUMENT .....  | 1         |
| 1.3. ACRONYMS AND ABBREVIATION .....  | 1         |
| 1.4. NOTE TO READER .....   | 2         |
| <b>2. IMPLEMENTATION MODEL</b> .....  | <b>3</b>  |
| 2.1. DAR-9400F / DAR-9500F .....  | 3         |
| 2.2. GATEWAY .....  | 4         |
| 2.3. APPLICATION DATA FLOW DIAGRAM .....  | 4         |
| 2.4. FUNCTIONAL DESCRIPTION OF AE'S .....   | 6         |
| 2.5. SEQUENCING OF REAL WORLD ACTIVITY .....  | 6         |
| <b>3. AE SPECIFICATIONS</b> .....   | <b>7</b>  |
| 3.1. <i>DAR-9400F/DAR-9500F</i> AE SPECIFICATION .....  | 7         |
| 3.2. <i>DAR-9400F/DAR-9500F</i> STORAGE MEDIA APPLICATION PROFILE CONFORMANCE STATEMENT ..... | 21        |
| <b>4. COMMUNICATION PROFILES</b> .....  | <b>22</b> |
| 4.1. TCP/IP STACK .....   | 22        |
| 4.2. EXTENSIONS/SPECIALIZATION/PRIVATIZATION .....  | 22        |
| 4.3. CONFIGURATION .....  | 22        |
| 4.4. SUPPORT FOR EXTENDED CHARACTER SETS .....  | 23        |
| <b>5. UID GENERATION</b> .....  | <b>24</b> |
| 5.1. TYPES OF GENERATED IMAGES .....  | 24        |
| 5.2. BASIC SPECIFICATION .....  | 25        |
| 5.3. RENUMBER SERIES BY SHOT ORDER MODE .....   | 27        |
| 5.4. UNIFIED SERIES NUMBER MODE .....   | 27        |
| 5.5. CHANGE SERIES INSTANCE UID MODE .....  | 28        |
| 5.6. SORT BY SHOT ORDER MODE .....  | 28        |
| <b>6. INFORMATION OBJECT DEFINITIONS</b> .....  | <b>29</b> |
| 6.1. X-RAY ANGIOGRAPHIC IMAGE .....   | 29        |
| 6.2. SECONDARY CAPTURE IMAGE .....  | 36        |
| 6.3. RDSR .....   | 40        |
| 6.4. DATA DICTIONARY OF PRIVATE ATTRIBUTES .....  | 50        |

This page is intentionally left blank.

## 1. Introduction

### 1.1. Purpose of This Document

The purpose of this document is to describe how *DAR-9400f/DAR-9500f* (with Gateway option) conforms to the DICOM standard. It describes what parts and definition it utilizes and in what way, in order to provide interoperability with other devices that claim same conformance.

### 1.2. Sources for This Document

American College of Cardiology –National Manufactures Association (ACR-NEMA) Digital Imaging and Communications V2.0  
ACR-NEMA Digital Imaging and Communications in Medicine (DICOM) v3.0, 2006.

### 1.3. Acronyms and Abbreviation

The following acronyms and abbreviations are used in this document.

- ACR American College of Radiology
- ACSE Association Control Service Element
- AE Application Entity
- ANSI American National Standards Institute
- AP Application Profile
- API Application Programming Interface
- ASCII American Standard Code for Information Interchange
- DICOM Digital Imaging and Communications in Medicine
- DIMSE DICOM Message Service Element
- DIMSE-C DICOM Message Service Element-Composite
- DIMSE-N DICOM Message Service Element-Normalized
- FSC File Set Creator
- FSR File Set Reader
- FSU File Set Updater
- GUI Graphical User Interface
- NEMA National Electrical Manufacturers Association
- PDU Protocol Data Unit
- RDSR DICOM Radiation Dose Structure Report
- RWA Real World Activity
- SCP Service Class Provider
- SCU Service Class User
- SOP Service Object Pair
- TCP/IP Transmission Control Protocol/Internet Protocol
- UID Unique Identifier
- MPPS Modality Performed Procedure Step
- PPS Performed Procedure Step
- PS Presentation State
- VM Value Multiplicity
- VR Value Representation
- VT Value Type

## 1.4. Note to Reader

- **Interoperability**

Interoperability refers to the ability of application functions, distributed over two or more systems, to work successfully together. The integration of medical devices into a networked environment may require application functions that are not specified within the scope of the DICOM standard. Consequently, using only the information provided by this conformance statement does not guarantee interoperability of Shimadzu Equipment with other vendor's equipment. It is the user's responsibility to thoroughly analyze the application requirements and to specify a solution that integrates Shimadzu equipment with the projected other vendor's equipment.

- **Validation**

Although Shimadzu equipment has been completely tested to verify that the implementation of the DICOM interface for this product corresponds with this Conformance Statement, even if comparison of respective Conformance Statement indicates that successful interconnection should be possible with another vendor's equipment, additional validation will always be necessary to ensure full functionality. It is the responsibility of the user to specify the appropriate test suite and to carry out the additional validation tests.

- **Version of the DICOM standard**

Shimadzu is committed to evolve with the DICOM standard as it adapts to meet the future requirement of users and technology. In order to do so, Shimadzu reserves the right to adapt and even discontinue delivery of its equipment. The user should ensure that any vendor whose equipment is connected to Shimadzu equipment also adapts to future version of the DICOM standard. If not, enhancement of Shimadzu may lead to loss of connectivity or interoperability.

- **Version Apply to *DAR-9400f/DAR-9500f***

Dicom Conformance Statement is applied to the following version of *DAR-9400f/DAR-9500f*. Refer to the old Dicom Conformance Statement to confirm older version of *DAR-9400f/DAR-9500f*.

- *DAR-9400f : Ver.4.2 or later*
- *DAR-9500f : Ver.6.0 or later*

## 2. Implementation Model

### 2.1. DAR-9400f / DAR-9500f

*DAR-9400f/DAR-9500f* is an acquisition and review station used for an angiographic image in the Cardiology environment.

The application, upon user request, will:

1. Acquire images from a CathLab and encapsulate them to the DICOM Standard Format.
2. Issue **C-ECHO** command to a remote DICOM SCP.
3. Issue **C-STORE** command to configured SCP in order to archive the acquired images.
4. Issue **C-MOVE** command to configured SCP.
5. Query (**C-FIND**), retrieve (**C-MOVE**) and display XA images from a remote DICOM SCP.
6. Query Modality Worklist (**C-FIND**)
7. Send **N-CREATE** and **N-SET** to MPPS server.
8. Read and display XA images from DICOM CD.
9. Act as FSC for DICOM CD. Write DICOM conformant CD-R
10. Act as FSR for DICOM CD. Read and display XA images from a DICOM conformant CD-R.
11. Print the images to remote DICOM SCP printer.

## 2.2. Gateway

Gateway is an option for *DAR-9400f/DAR-9500f* image transfer system. *DAR-9400f/DAR-9500f* sends the image to the Gateway, edit the image properly and send it to the set SCP.

## 2.3. Application Data Flow Diagram

### 2.3.1. DAR-9400f / DAR-9500f

(1) Verification

*DAR-9400f/DAR-9500f* can send **C-ECHO** verifications.

(2) Basic Worklist Management or IHE Worklist Management

*DAR-9400f/DAR-9500f* will issue a **C-FIND** for the hospital Worklist.

(3) Find

*DAR-9400f/DAR-9500f* will issue a **C-FIND** command to a remote SCP to retrieve information about the studies stored on the remote SCP.

(4) Move Images

*DAR-9400f/DAR-9500f* will issue a **C-MOVE** command to a remote SCP to copy study information from one SCP to another or from a remote SCP to itself.

(5) Store Images and RDSR as SCP

*DAR-9400f/DAR-9500f* will receive process and accept **C-STORE** command from a remote SCU and if the association succeeds, it will store the received data on its physical storage space.

(6) Store Images and RDSR as SCU

*DAR-9400f/DAR-9500f* or Gateway will issue a **C-STORE** command to a remote SCP. If the association is successful it will send images for storage on the remote SCP.

(7) MPPS Management

*DAR-9400f/DAR-9500f* will issue an **N-CREATE** event to notify the creation of a new acquisition study and issue an **N-SET** event when this acquisition study is completed.

(8) Print

*DAR-9400f/DAR-9500f* will print an image or loop of images to the remote DICOM SCP printer.


(9) Read CD-R

*DAR-9400f/DAR-9500f* will read any DICOM conformant CD-R although it will only display compatible images.

(10) Write CD-R

*DAR-9400f/DAR-9500f* will write a DICOM conformant CD-R for the supported SOP classes.





## 2.4. Functional Description of AE's

### (1) DAR-9400f / DAR-9500f

The *DAR-9400f/DAR-9500f* AE acts as a SCU and a SCP.

- If configured, *DAR-9400f/DAR-9500f* can query for the patient Worklist. The list of scheduled patient will be presented to the user and all fields in the patient demographic entry forms will be filled with the chosen patient. If all mandatory fields cannot be filled a form will be presented to the user with the missing fields highlighted. *DAR-9400f/DAR-9500f* will issue a **C-Find** request to retrieve Worklist information from a remote Modality Worklist SCP.
- *DAR-9400f/DAR-9500f* is a system designed to acquire images coming from the FPD of a catheterization laboratory. The system then compresses these images and transmits these compressed images over the network to a remote server that will place them in a safe place for long-term archive and retrieval.
- If configured to do so, when a study is terminated, *DAR-9400f/DAR-9500f* will send a command to the remote server to move images and RDSR to one or more secondary destinations. The move may involve all the images objects of a study or only specific images belonging to one or more series. There may be any number of secondary destinations. When "Secondary Destination" is/are configured, *DAR-9400f/DAR-9500f* will issue a C-MOVE command to the server configured as "Primary Server" using either "Move by study" or "Move by Series" UID. The command will contain the "Primary Server" as origin and the "Secondary destination" as destination.
- *DAR-9400f/DAR-9500f* will issue a C-FIND command to get and display the content of a STORAGE-SCP remote server; it will then issue a C-MOVE command at study root level to the same server to transfer the selected image data proposing itself as destination.
- *DAR-9400f/DAR-9500f* will issue DIMSE-N set command to print SCP to print the image.
- *DAR-9400f/DAR-9500f* will issue a C-STORE command to a configured remote SCP to store images previously read from a DICOM CD.
- *DAR-9400f/DAR-9500f* will read data from a DIOCM CD that is present in its CD drive when requested to do so by the user
- *DAR-9400f/DAR-9500f* will accept association from remote SCU and accept and process C-STORE commands for DICOM Data Object of the allowed SOP classes.
- *DAR-9400f/DAR-9500f* will issue N-CREATE and N-SET command to MPPS server to notice the implementation of study.

## 2.5. Sequencing of Real World Activity

The storage Verification is done when a study is closed, and only if files for the study have been transmitted for storage to a remote SCP.

The Storage verification is done after the current study is closed.

The physical CD-R writing can only occur after an empty CD-R is inserted in the drive.

### 3. AE Specifications

#### 3.1. DAR-9400f/DAR-9500f AE Specification

*DAR-9400f/DAR-9500f* provides Standard Conformance to the following DICOM V3.0 SOP Class as an SCU.

**Table 1 Verification SOP Class as SCU**

| SOP Class Name                        | SOP Class UID                 |
|---------------------------------------|-------------------------------|
| Verification                          | 1.2.840.10008.1.1             |
| Study Root Query/Retrieve IM Find     | 1.2.840.10008.5.1.4.1.2.2.1   |
| Study Root Query/Retrieve IM Move     | 1.2.840.10008.5.1.4.1.2.2.2   |
| XA – X-ray Angiographic image storage | 1.2.840.10008.5.1.4.1.1.12.1  |
| Secondary Capture Image storage       | 1.2.840.10008.5.1.4.1.1.7     |
| Modality Worklist SOP class           | 1.2.840.10008.5.1.4.31        |
| Basic Grayscale Print Management Meta | 1.2.840.10008.5.1.1.9         |
| X-Ray Radiation Dose SR               | 1.2.840.10008.5.1.4.1.1.88.67 |

### 3.1.1. Association Establishment Policies

#### 3.1.1.1. General

The following Application Context Name will be proposed and recognized by *DAR-9400f/DAR-9500f*.

- DICOM 3.0 Application Context      **1.2.840.10008.3.1.1.1**

#### 3.1.1.2. Number of Associations

The maximum number of association accepted or maintained by *DAR-9400f/DAR-9500f* is limited only by the physical memory of the machine on which it runs. Typically it can be up to 10. However, only 1 of it will be sent from Gateway.

#### 3.1.1.3. Asynchronous Nature

*DAR-9400f/DAR-9500f* allows a single outstanding operation on any association. Therefore, *DAR-9400f/DAR-9500f* does not support asynchronous operations window negotiation, other than the default as specified by the DICOM specification.

#### 3.1.1.4. Implementation Identifying Information

*DAR-9400f/DAR-9500f* will respond with the following implementation identifying parameters:

- For Safire  
Implementation Class UID (acquisition) **1.2.392.200036.9110.13.XXXXXXXXXXXXXXX**  
Implementation Class UID (review)      **1.2.392.200036.9110.14.XXXXXXXXXXXXXXX**
- For Alexa  
Implementation Class UID (acquisition) **1.2.392.200036.9110.15.XXXXXXXXXXXXXXX**  
Implementation Class UID (review)      **1.2.392.200036.9110.16.XXXXXXXXXXXXXXX**
- For Trinias  
Implementation Class UID (acquisition) **1.2.392.200036.9110.17.XXXXXXXXXXXXXXX**  
Implementation Class UID (review)      **1.2.392.200036.9110.18.XXXXXXXXXXXXXXX**

The last number of the implementation class UID is the 13 digits maximum machine serial number.

Implementation Version Name      **VOYAGER\_VX\_X\_X**  
Where X\_X\_X is the software version

The implementation version name policies are the following: product name “**VOYAGER**” followed by the version of the product, “**\_v1\_0\_0**”.

When send the image via Gateway, it will respond with the following implementation identification parameters: 1.2.392.200036.9110.1.0.6711.2000731

Implementation version name is: SRVP-REC2.10.

### 3.1.2. Association Initiation by Real World Activity

#### 3.1.2.1. Real World Activity – Verification

(1) Associated Real World Activity – Verification

*DAR-9400f/DAR-9500f* will send C-ECHO to verify the other systems if they are activated.

(2) Presentation context Table – Verification

*DAR-9400f/DAR-9500f* supports the transfer syntaxes listed in Table 2. For a **Verification** request, *DAR-9400f/DAR-9500f* supports the Presentation Contexts listed in Table 2.

**Table 2 Presentation Contexts**

| Abstract Syntax |                   | Transfer Syntax              |                   | Role | Extended Negotiation |
|-----------------|-------------------|------------------------------|-------------------|------|----------------------|
| SOP Class       | SOP Class UID     | Name List                    | UID List          |      |                      |
| Verification    | 1.2.840.10008.1.1 | Implicit VR<br>Little Endian | 1.2.840.10008.1.2 | SCU  | None                 |

(3) SOP Specific conformance – Verification

*DAR-9400f/DAR-9500f* provides the standard conformance to the DICOM Verification SOP class.

#### 3.1.2.2. Real World Activity – Query Study

(1) Associated Real World Activity – Query Study

*DAR-9400f/DAR-9500f* will issue a **FIND** request when a user of *DAR-9400f/DAR-9500f/DAR-9500f* wishes to view patient and study information from a remote DICOM SCP.

(2) Presentation context Table – Query Study

*DAR-9400f/DAR-9500f* supports the transfer syntaxes listed in Table 3. For a **QUERY** request, *DAR-9400f/DAR-9500f* supports the Presentation Contexts listed in Table 3.

**Table 3 Presentation Contexts**

| Abstract Syntax   |                             | Transfer Syntax              |                   | Role | Extended Negotiation |
|---|-----------------------------|------------------------------|-------------------|------|----------------------|
| SOP Class   | SOP Class UID               | Name List                    | UID List          |      |                      |
| Study Root<br>Query/Retrieve<br>Information Model<br>– FIND | 1.2.840.10008.5.1.4.1.2.2.1 | Implicit VR<br>Little Endian | 1.2.840.10008.1.2 | SCU  | None                 |

(3) SOP Specific conformance – Query Study

*DAR-9400f/DAR-9500f* uses Relational query with Study root level by default. If the extended negotiation is not successful, it uses Relational query with Patient root level by default. If the extended negotiation still does not succeed, *DAR-9400f/DAR-9500f* will use Hierarchical query with patient root model.

➤ Matching key attribution

*DAR-9400f/DAR-9500f* Request matching of the following key attributes.

**Table 4 Matching Key Attributes - Study Root Query/Retrieve Information Model \***

| Attribute Name      | Type | Tag         |
|---------------------|------|-------------|
| <b>STUDY LEVEL</b>  |      |             |
| Study Date          | R    | (0008,0020) |
| Accession Number    | R    | (0008,0050) |
| Patient's Name      | R    | (0010,0010) |
| Patient ID          | R    | (0010,0020) |
| <b>SERIES LEVEL</b> |      |             |
| Modality            | R    | (0008,0060) |

\* User can input these items on the monitor as a search key.

**3.1.2.3. Real World Activity - Move Images**

(1) Associated Real World Activity – Move Images

*DAR-9400f/DAR-9500f* will issue a MOVE request when a user of *DAR-9400f/DAR-9500f* wishes to move one or more studies from a remote DICOM SCP back to *DAR-9400f/DAR-9500f* (retrieve) or another remote DICOM SCP.

(2) Presentation context Table – Move

*DAR-9400f/DAR-9500f* supports the transfer syntaxes listed in Table 5. For a **MOVE** request, *DAR-9400f/DAR-9500f* supports the Presentation Contexts listed in Table 5 and Table 6.

**Table 5 Move Transfer Syntaxes**

| Transfer Syntax                                 | UID               |
|---|-------------------|
| DICOM Implicit VR Little Endian Transfer Syntax | 1.2.840.10008.1.2 |

**Table 6 Move Presentation Contexts**

| Abstract Syntax                   |                             | Transfer Syntax                 |                   | Role | Extended Negotiation |
|-----------------------------------|-----------------------------|---------------------------------|-------------------|------|----------------------|
| SOP Class                         | SOP Class UID               | Name List                       | UID List          |      |                      |
| Study Root Query/Retrieve IM Move | 1.2.840.10008.5.1.4.1.2.2.2 | DICOM Implicit VR Little Endian | 1.2.840.10008.1.2 | SCU  | None                 |

(3) SOP Specific Conformance – Move

*DAR-9400f/DAR-9500f* uses specific keys for Move operation. When doing a series move the Study UID and Series UID are used as keys. When doing a study move only the Study UID is used as key.

### 3.1.2.4. Real World Activity – Request to Store Images and Dose Information

#### (1) Associated Real World Activity – Storage as SCU

*DAR-9400f/DAR-9500f* will issue a **Storage** request when a user of *DAR-9400f/DAR-9500f* wishes to send a study of images to a remote DICOM SCP.

#### (2) Presentation context Table – Storage as SCU

*DAR-9400f/DAR-9500f* supports the transfer syntaxes listed in Table 7.

**Table 7 Worklist Presentation Context**

| Abstract Syntax                        |                               | Transfer Syntax   |                        | Role | Extended Negotiation |
|--|-------------------------------|---|------------------------|------|----------------------|
| SOP Class                              | SOP Class UID                 | Name List   | UID List               |      |                      |
| X-ray<br>Angiographic<br>Image Storage | 1.2.840.10008.5.1.4.1.1.12.1  | Implicit VR<br>Little Endian                              | 1.2.840.10008.1.2      | SCU  | None                 |
|  |                               | Explicit VR<br>Little Endian                              | 1.2.840.10008.1.2.1    |      |                      |
|  |                               | JPEG Lossless,<br>hierarchical, first<br>order prediction | 1.2.840.10008.1.2.4.70 |      |                      |
| Secondary<br>Capture Image<br>Storage  | 1.2.840.10008.5.1.4.1.1.7     | Implicit VR<br>Little Endian                              | 1.2.840.10008.1.2      | SCU  | None                 |
|  |                               | Explicit VR<br>Little Endian                              | 1.2.840.10008.1.2.1    |      |                      |
|  |                               | JPEG Lossless,<br>hierarchical, first<br>order prediction | 1.2.840.10008.1.2.4.70 |      |                      |
| X-Ray<br>Radiation<br>Dose SR          | 1.2.840.10008.5.1.4.1.1.88.67 | Implicit VR Little<br>Endian                              | 1.2.840.10008.1.2      | SCU  | None                 |

#### (3) Storage Presentation Contexts selection

- Transfer syntax can be configured. They can be enabled/disabled and the presented order can be selected
- If no Transfer syntax are selected in the option then the system try to negotiate the default Transfer Syntax (Original Storage Transfer Syntaxes) and Implicit VR Little Endian Transfer Syntax.

### 3.1.2.5. Real World Activity – Query Worklist

(1) Associated Real World Activity – Query Worklist as SCU

*DAR-9400f/DAR-9500f* will issue a **query Worklist** request when a user of *DAR-9400f/DAR-9500f* opens a new study if a Modality Worklist SCP is configured in its host table.

(2) Presentation context Table – Query Worklist as SCU

*DAR-9400f/DAR-9500f* supports the transfer syntaxes listed in Table 8. For a **Query Worklist** request, *DAR-9400f/DAR-9500f* supports the Presentation Contexts listed in Table 8.

**Table 8 Worklist Presentation Contexts**

| Abstract Syntax                            |                        | Transfer Syntax                 |                   | Role | Extended Negotiation |
|--|------------------------|---------------------------------|-------------------|------|----------------------|
| SOP Class                                  | SOP Class UID          | Name List                       | UID List          |      |                      |
| Modality Worklist Information Model - FIND | 1.2.840.10008.5.1.4.31 | DICOM Implicit VR Little Endian | 1.2.840.10008.1.2 | SCU  | None                 |

(3) SOP Specific Conformance – general purpose Worklist SOP Class as SCU

*DAR-9400f/DAR-9500f* supports queries against the Worklist Information Model using the baseline **C-FIND SCU** behaviour.

Also, *DAR-9400f/DAR-9500f* supports the character sets that are indicated on section “4.4Support for Extended Character Sets”.

- *DAR-9400f/DAR-9500f* Request matching of the following key attributes

**Table 9 Matching Key Attributes - Query Worklist\***

| Attribute Name                         | VR | Tag          |
|--|----|--------------|
| <b>Scheduled Procedure Step Module</b> |    |              |
| Scheduled Procedure Step Sequence      | SQ | (0040,0100)  |
| >Modality**                            | CS | >(0008,0060) |
| >Scheduled Station AE Title            | AE | >(0040,0001) |
| >Scheduled Procedure Step Start Date   | DA | >(0040,0002) |
| >Scheduled Performing Physician's Name | PN | >(0040,0006) |
| <b>Requested Procedure Module</b>      |    |              |
| Requested Procedure ID                 | SH | (0040,1001)  |
| <b>Imaging Service Request Module</b>  |    |              |
| Accession Number                       | SH | (0008,0050)  |
| <b>Patient Identification Module</b>   |    |              |
| Patient's Name                         | PN | (0010,0010)  |
| Patient ID                             | LO | (0010,0020)  |

\* User can input these items on the monitor as a search key.

\*\* Always searching with XA.

- *DAR-9400f/DAR-9500f* supports the character sets listed in Table 24.



### 3.1.2.6. Real World Activity – Request to Print Images

(1) Associated Real World Activity

*DAR-9400f/DAR-9500f* will issue a **Print** request when the user wants to send study images to the remote DICOM printer SCP.

(2) Presentation context Table – Request to Print Images

*DAR-9400f/DAR-9500f* supports the transfer syntaxes listed in Table 10. For a **Print** request, *DAR-9400f/DAR-9500f* supports the Presentation Contexts listed in Table 10.

**Table 10 Request to Print Images Presentation Contexts**

| Abstract Syntax                       |                       | Transfer Syntax           |                     | Role | Extended Negotiation |
|---------------------------------------|-----------------------|---------------------------|---------------------|------|----------------------|
| SOP Class                             | SOP Class UID         | Name List                 | UID List            |      |                      |
| Basic Grayscale Print Management Meta | 1.2.840.10008.5.1.1.9 | Explicit VR Little Endian | 1.2.840.10008.1.2.1 | SCU  | None                 |

(3) SOP Specific Conformance

*DAR-9400f/DAR-9500f* provides the standard conformance to the DICOM Basic Grayscale Print Management Meta SOP class.

### 3.1.2.7. Real World Activity – Create a new acquisition study

(1) Associated Real World Activity – create a new acquisition study

*DAR-9400f/DAR-9500f* will issue an MPPS N-CREATE event when a user of *DAR-9400f/DAR-9500f* creates a new study in acquisition, if a PPS Manager is configured in its host table.

(2) Presentation context Table – MPPS N-CREATE

*DAR-9400f/DAR-9500f* supports the transfer syntaxes listed in Table 11 for MPPS N-CREATE event, *DAR-9400f/DAR-9500f* supports the Presentation Contexts listed in Table 12.

**Table 11 MPPS N-CREATE Transfer Syntaxes**

| Transfer Syntax                                 | UID               |
|---|-------------------|
| DICOM Implicit VR Little Endian Transfer Syntax | 1.2.840.10008.1.2 |

**Table 12 MPPS N-CREATE Presentation Contexts**

| Abstract Syntax                             |                         | Transfer Syntax   | Role | Extended Negotiation |
|---|-------------------------|-------------------|------|----------------------|
| SOP Class                                   | SOP Class UID           |                   |      |                      |
| Modality Performed Procedure Step SOP class | 1.2.840.10008.3.1.2.3.3 | all from Table 11 | SCU  | None                 |

(3) SOP Specific Conformance

DAR-9400f/DAR-9500f provides the following table describes the supported attributes of a N-CREATE message.

**Table 13 Performed Procedure Step N-CREATE Attributes**

| Tag   | VR | Attribute Name  | Value of N-CREATE                            |
|---|----|---|--|
| <b>SOP Common Module</b>                            |    |   |  |
| (0008,0005)   | CS | Specific Character Set  | "ISO IR 6"<br>"ISO IR100"<br>"ISO 2022 IR87" |
| <b>Performed Procedure Step Relationship Module</b> |    |   |  |
| (0010,0010)   | PN | Patient's Name  | From Modality Worklist or user input         |
| (0010,0020)   | LO | Patient ID  | From Modality Worklist or user input         |
| (0010,0030)   | DA | Patient's Birth Date  | From Modality Worklist or user input         |
| (0010,0040)   | CS | Patient's Sex   | From Modality Worklist or user input         |
| (0008,1120)   | SQ | Referenced Patient Sequence                                   | Zero length                                  |
| >(0008,1150)  | UI | Referenced SOP Class UID                                      | N.A.   |
| >(0008,1155)  | UI | Referenced SOP Instance UID                                   | N.A.   |
| >(0020,000D)  | UI | Study Instance UID  | From Modality Worklist                       |
| >(0008,1110)  | SQ | Referenced Study Sequence                                     | Zero length                                  |
| >(0008,0050)  | SH | Accession Number  | From Modality Worklist or user input         |
| >(0032,1060)  | LO | Requested Procedure Description                               | From Modality Worklist                       |
| >(0040,0009)  | SH | Scheduled Procedure Step ID                                   | From Modality Worklist                       |
| >(0040,0007)  | LO | Scheduled Procedure Step Description                          | From Modality Worklist                       |
| >(0040,0008)  | SQ | Scheduled Protocol Code Sequence                              | Zero length                                  |
| (0040,0270)   | SQ | Scheduled Step Attributes Sequence                            | Zero length                                  |
| >(0040,1001)  | SH | Requested Procedure ID  | From Modality Worklist                       |
| <b>Performed Procedure Step Information Module</b>  |    |   |  |
| (0040,0241)   | AE | Performed Station AE Title                                    | MPPS AE Title                                |
| (0040,0242)   | SH | Performed Station Name  | From configuration                           |
| (0040,0243)   | SH | Performed Location  | Zero length                                  |
| (0040,0244)   | DA | Performed Procedure Step Start Time                           | Actual start date                            |
| (0040,0245)   | TM | Performed Procedure Step Start Time                           | Actual start time                            |
| (0040,0250)   | DA | Performed Procedure Step End Date                             | Zero length                                  |
| (0040,0251)   | TM | Performed Procedure Step End Time                             | Zero length                                  |
| (0040,0252)   | CS | Performed Procedure Step Status                               | IN PROGRESS                                  |
| (0040,0253)   | SH | Performed Procedure Step ID                                   | Automatically created                        |
| (0040,0254)   | LO | Performed Procedure Step Description                          | Zero length                                  |
| (0040,0255)   | LO | Performed Procedure Type Description                          | Zero length                                  |
| (0008,1032)   | SQ | Procedure Code Sequence                                       | Zero length                                  |
| (0040,0281)   | SQ | Performed Procedure Step Discontinuation Reason Code Sequence | N.A.   |
| >(0008,0100)  | SH | Code Value  | N.A.   |
| >(0008,0102)  | SH | Coding Scheme Designator                                      | N.A.   |
| >(0008,0104)  | LO | Code Meaning  | N.A.   |

| Tag                                    | VR | Attribute Name                               | Value of N-CREATE  |
|--|----|--|--|
| <b>Image Acquisition Result Module</b> |    |  |  |
| (0008,0060)                            | CS | Modality                                     | XA   |
| (0020,0010)                            | SH | Study ID                                     | From Modality Worklist or user input                           |
| (0040,0260)                            | SQ | Performed Protocol Code Sequence             | Zero or more items   |
| >(0008,0100)                           | SH | Code Value                                   | An appropriate value will be sent if Sequence Item is present. |
| >(0008,0102)                           | SH | Coding Scheme Designator                     | An appropriate value will be sent if Sequence Item is present. |
| >(0008,0103)                           | SH | Coding Scheme Version                        | An appropriate value will be sent if Sequence Item is present. |
| >(0008,0104)                           | LO | Code Meaning                                 | An appropriate value will be sent if Sequence Item is present. |
| (0040,0340)                            | SQ | Performed Series Sequence                    | Zero length  |
| >(0008,1050)                           | PN | Performing Physician's Name                  | N.A.   |
| >(0008,1070)                           | PN | Operator's Name                              | N.A.   |
| >(0018,1030)                           | LO | Protocol Name                                | N.A.   |
| >(0020,000E)                           | UI | Series Instance UID                          | N.A.   |
| >(0008,103E)                           | LO | Series Description                           | N.A.   |
| >(0008,0054)                           | AE | Retrieve AE Title                            | N.A.   |
| >(0008,1140)                           | SQ | Referenced Image Sequence                    | N.A.   |
| >(0040,0220)                           | SQ | Referenced Standalone SOP Instance Seq.      | N.A.   |
| <b>Radiation Dose Module</b>           |    |  |  |
| (0008,2229)                            | SQ | Anatomic Structure, Space or Region Sequence | No length  |
| >(0008,0100)                           | SH | Code Value                                   | N.A.   |
| >(0008,0102)                           | SH | Coding Scheme Designator                     | N.A.   |
| >(0008,0104)                           | LO | Code Meaning                                 | N.A.   |
| (0040,0300)                            | US | Total Time of Fluoroscopy                    | Zero length  |
| (0040,0301)                            | US | Total Number of Exposures                    | Zero length  |
| (0040,0302)                            | US | Entrance Dose                                | Zero length  |
| (0040,8302)                            | DS | Entrance Dose in mGy                         | Zero length  |
| (0018,115E)                            | DS | Image and Fluoroscopy Area Dose Product      | Zero length  |
| (0040,0310)                            | ST | Comments on Radiation Dose                   | Zero length  |
| (0040,030E)                            | SQ | Exposure Dose Sequence                       | Zero length  |
| >(0018,115A)                           | CS | Radiation Mode                               | N.A.   |
| >(0018,0060)                           | DS | KVp  | N.A.   |
| >(0018,8151)                           | DS | X-Ray Tube Current in $\mu$ A                | N.A.   |
| >(0018,1150)                           | IS | Exposure Time                                | N.A.   |
| >(0018,1160)                           | SH | Filter Type                                  | N.A.   |
| >(0018,7050)                           | CS | Filter Material                              | N.A.   |
| >(0040,0310)                           | ST | Comments on Radiation Dose                   | N.A.   |

Enable to send the following attributes from *DAR-9400f* Ver.4.2 or later version and *DAR-9500f* Ver.5.2.40 or later version.

### 3.1.2.8. Real World Activity – Close a study in acquisition

(1) Associated Real World Activity – close a study in acquisition

*DAR-9400f/DAR-9500f* will issue an MPPS N-SET event when a user of *DAR-9400f/DAR-9500f* closes a study in acquisition, if a PPS Manager is configured in its host table.

(2) Presentation context Table – MPPS N-SET

*DAR-9400f/DAR-9500f* supports the transfer syntaxes listed in Table 14 for MPPS N-SET event;  
*DAR-9400f/DAR-9500f* supports the Presentation Contexts listed in Table 15.

**Table 14 MPPS N-SET Transfer Syntaxes**

| Transfer Syntax                                 | UID               |
|---|-------------------|
| DICOM Implicit VR Little Endian Transfer Syntax | 1.2.840.10008.1.2 |

**Table 15 MPPS N-CREATE Presentation Contexts**

| Abstract Syntax                             |                         | Transfer Syntax   | Role | Extended Negotiation |
|---|-------------------------|-------------------|------|----------------------|
| SOP Class                                   | SOP Class UID           |                   |      |                      |
| Modality Performed Procedure Step SOP class | 1.2.840.10008.3.1.2.3.3 | all from Table 14 | SCU  | None                 |

(3) SOP Specific Conformance

DAR-9400f/DAR-9500f provides The following table describes the supported attributes of a N-SET message.

**Table 16 Performed Procedure Step N-SET Attributes**

| Tag   | VR | Attribute Name  | Value of N-SET   |
|---|----|---|--|
| <b>SOP Common Module</b>                            |    |   |  |
| (0008,0005)   | CS | Specific Character Set  | N.A.   |
| <b>Performed Procedure Step Relationship Module</b> |    |   |  |
| (0010,0010)   | PN | Patient's Name  | N.A.   |
| (0010,0020)   | LO | Patient ID  | N.A.   |
| (0010,0030)   | DA | Patient's Birth Date  | N.A.   |
| (0010,0040)   | CS | Patient's Sex   | N.A.   |
| (0008,1120)   | SQ | Referenced Patient Sequence                                   | N.A.   |
| >(0008,1150)  | UI | Referenced SOP Class UID                                      | An appropriate value will be sent                              |
| >(0008,1155)  | UI | Referenced SOP Instance UID                                   | An appropriate value will be sent                              |
| >(0020,000D)  | UI | Study Instance UID  | N.A.   |
| >(0008,1110)  | SQ | Referenced Study Sequence                                     | N.A.   |
| >(0008,0050)  | SH | Accession Number  | N.A.   |
| >(0032,1060)  | LO | Requested Procedure Description                               | N.A.   |
| >(0040,0009)  | SH | Scheduled Procedure Step ID                                   | N.A.   |
| >(0040,0007)  | LO | Scheduled Procedure Step Description                          | N.A.   |
| >(0040,0008)  | SQ | Scheduled Protocol Code Sequence                              | N.A.   |
| (0040,0270)   | SQ | Scheduled Step Attributes Sequence                            | N.A.   |
| >(0040,1001)  | SH | Requested Procedure ID  | N.A.   |
| <b>Performed Procedure Step Information Module</b>  |    |   |  |
| (0040,0241)   | AE | Performed Station AE Title                                    | N.A.   |
| (0040,0242)   | SH | Performed Station Name  | N.A.   |
| (0040,0243)   | SH | Performed Location  | N.A.   |
| (0040,0244)   | DA | Performed Procedure Step Start Time                           | N.A.   |
| (0040,0245)   | TM | Performed Procedure Step Start Time                           | N.A.   |
| (0040,0250)   | DA | Performed Procedure Step End Date                             | Actual end date  |
| (0040,0251)   | TM | Performed Procedure Step End Time                             | Actual end time  |
| (0040,0252)   | CS | Performed Procedure Step Status                               | DISCONTINUED or COMPLETED                                      |
| (0040,0253)   | SH | Performed Procedure Step ID                                   | N.A.   |
| (0040,0254)   | LO | Performed Procedure Step Description                          | N.A.   |
| (0040,0255)   | LO | Performed Procedure Type Description                          | N.A.   |
| (0008,1032)   | SQ | Procedure Code Sequence                                       | N.A.   |
| (0040,0281)   | SQ | Performed Procedure Step Discontinuation Reason Code Sequence | Zero length  |
| >(0008,0100)  | SH | Code Value  | An appropriate value will be sent if Sequence Item is present. |
| >(0008,0102)  | SH | Coding Scheme Designator                                      | An appropriate value will be sent if Sequence Item is present. |
| >(0008,0104)  | LO | Code Meaning  | An appropriate value will be sent if Sequence Item is present. |

| Tag                                    | VR | Attribute Name                               | Value of N-SET   |
|--|----|--|--|
| <b>Image Acquisition Result Module</b> |    |  |  |
| (0008,0060)                            | CS | Modality                                     | N.A.   |
| (0020,0010)                            | SH | Study ID                                     | N.A.   |
| (0040,0260)                            | SQ | Performed Protocol Code Sequence             | Zero or more items   |
| >(0008,0100)                           | SH | Code Value                                   | An appropriate value will be sent if Sequence Item is present. |
| >(0008,0102)                           | SH | Coding Scheme Designator                     | An appropriate value will be sent if Sequence Item is present. |
| >(0008,0103)                           | SH | Coding Scheme Version                        | An appropriate value will be sent if Sequence Item is present. |
| >(0008,0104)                           | LO | Code Meaning                                 | An appropriate value will be sent if Sequence Item is present. |
| (0040,0340)                            | SQ | Performed Series Sequence                    | One or more items  |
| >(0008,1050)                           | PN | Performing Physician's Name                  | An appropriate value will be sent                              |
| >(0008,1070)                           | PN | Operator's Name                              | Zero length  |
| >(0018,1030)                           | LO | Protocol Name                                | An appropriate value will be sent                              |
| >(0020,000E)                           | UI | Series Instance UID                          | An appropriate value will be sent                              |
| >(0008,103E)                           | LO | Series Description                           | An appropriate value will be sent                              |
| >(0008,0054)                           | AE | Retrieve AE Title                            | An appropriate value will be sent                              |
| >(0008,1140)                           | SQ | Referenced Image Sequence                    | One or more items.   |
| >(0040,0220)                           | SQ | Referenced Standalone SOP Instance Seq.      | Zero length  |
| <b>Radiation Dose Module</b>           |    |  |  |
| (0008,2229)                            | SQ | Anatomic Structure, Space or Region Sequence | One or more items  |
| >(0008,0100)                           | SH | Code Value                                   | An appropriate value will be sent if Sequence Item is present. |
| >(0008,0102)                           | SH | Coding Scheme Designator                     | An appropriate value will be sent if Sequence Item is present. |
| >(0008,0104)                           | LO | Code Meaning                                 | An appropriate value will be sent if Sequence Item is present. |
| (0040,0300)                            | US | Total Time of Fluoroscopy                    | Actual total time of fluoroscopy                               |
| (0040,0301)                            | US | Total Number of Exposures                    | Actual total number of exposures                               |
| (0040,0302)                            | US | Entrance Dose                                | Actual entrance Dose in dGy                                    |
| (0040,8302)                            | DS | Entrance Dose in mGy                         | Actual entrance Dose in mGy                                    |
| (0018,115E)                            | DS | Image and Fluoroscopy Area Dose Product      | Actual image area dose product                                 |
| (0040,0310)                            | ST | Comments on Radiation Dose                   |  |
| (0040,030E)                            | SQ | Exposure Dose Sequence                       | One or more items  |
| >(0018,115A)                           | CS | Radiation Mode                               | Specified X-Ray radiation mode                                 |
| >(0018,0060)                           | DS | KVp  | Peak kilo voltage output of the x-ray generator                |
| >(0018,8151)                           | DS | X-Ray Tube Current in $\mu$ A                | X-Ray tube current in $\mu$ A                                  |
| >(0018,1150)                           | IS | Exposure Time                                | The time of X-Ray exposure                                     |
| >(0018,1160)                           | SH | Filter Type                                  | Type of filter(s)  |
| >(0018,7050)                           | CS | Filter Material                              | The X-Ray absorbing material used in the filter                |
| >(0040,0310)                           | ST | Comments on Radiation Dose                   | User-defined comments on any special conditions                |

Enable to send the following attributes from *DAR-9400f* Ver.4.2 or later version and *DAR-9500f* Ver.5.2.40 or later version.

### 3.1.3. Association Acceptance Policy

#### 3.1.3.1. Real World Activity - Store Received Images and Dose Information

(1) Associated Real World Activity – Store Received Images

*DAR-9400f/DAR-9500f* will archive images that are sent from **C-STORE SCU**.

(2) Presentation Context Table – Store Received Images

*DAR-9400f/DAR-9500f* supports the following transfer syntaxes listed in Table 17.

**Table 17 Storage Transfer Syntaxes**

| Abstract Syntax                        |                               | Transfer Syntax  |                        | Role | Extended Negotiation |
|--|-------------------------------|--|------------------------|------|----------------------|
| SOP Class                              | SOP Class UID                 | Name List  | UID List               |      |                      |
| X-ray<br>Angiographic<br>Image Storage | 1.2.840.10008.5.1.4.1.1.12.1  | Implicit VR<br>Little Endian                                 | 1.2.840.10008.1.2      | SCP  | None                 |
|  |                               | Explicit VR<br>Little Endian                                 | 1.2.840.10008.1.2.1    |      |                      |
|  |                               | JPEG Lossless,<br>hierarchical,<br>first<br>order prediction | 1.2.840.10008.1.2.4.70 |      |                      |
| Secondary<br>Capture Image<br>Storage  | 1.2.840.10008.5.1.4.1.1.7     | Implicit VR<br>Little Endian                                 | 1.2.840.10008.1.2      | SCP  | None                 |
|  |                               | Explicit VR<br>Little Endian                                 | 1.2.840.10008.1.2.1    |      |                      |
|  |                               | JPEG Lossless,<br>hierarchical,<br>first<br>order prediction | 1.2.840.10008.1.2.4.70 |      |                      |
| X-Ray<br>Radiation<br>Dose SR          | 1.2.840.10008.5.1.4.1.1.88.67 | Implicit VR<br>Little Endian                                 | 1.2.840.10008.1.2      | SCP  | None                 |

*DAR-9400f/DAR-9500f* returns one of the following status codes (Table 18).

**Table 18 Storage status codes**

| Service Status | Further Meaning                   | Protocol Codes | Related Fields | Description  |
|----------------|-----------------------------------|----------------|----------------|--|
| Refused        | Out of resources                  | A700           |                | Indicates that there was not enough storage space to store the image. Recovery from this condition is left to the administrative functions available in <i>DAR-9400f/DAR-9500f/DAR-9500f</i> . |
|                | SOP Class not supported           | A800           |                | Indicates that the SOP Class of the Image in the <b>C-Store</b> operation did not match the Abstract Syntax negotiated for the Presentation Context.   |
| Error          | Data set does not match SOP Class | A900           |                | Indicates that the Data Set does not encode an instance of the SOP Class specified.  |
|                | Failed                            | C000           |                | The operation was not successful.  |
|                | Cannot understand                 | C005           |                | Indicates that the Data Set cannot be parsed into elements by <i>DAR-9400f/DAR-9500f/DAR-9500f</i> .   |
| Warning        | Coercion of data elements         | B000           |                | Data elements were modified before being stored.   |
|                | Data set does not match SOP Class | B007           |                | Indicates that the Data Set does not match the SOP Class, but that the image was stored anyway.  |
|                | Elements Discarded                | B006           |                | Indicates that some of the elements of the Data Set were discarded.  |
|                | Duplicate SOP Instance UID        | D000           |                | Indicates that the SOP Instance UID of the specified image is already stored in the database.  |
| Success        | Success                           | 0000           |                | Operation performed properly.  |

### (3) Presentation Context Acceptance Criterion – Store Received Images

*DAR-9400f/DAR-9500f* will accept any number of **Storage** Presentation Contexts per association request. Any one Abstract Syntax may be specified more than once in an association request, if the Transfer Syntaxes differ between the Presentation Contexts.

The acceptable Presentation Contexts which *DAR-9400f/DAR-9500f* may accept are specified in Table 15. *DAR-9400f/DAR-9500f* will examine proposed Presentation Contexts in the order proposed. The first acceptable Presentation Context (other than Verification) determines the Abstract Syntax which will be used for the association.



### 3.2. DAR-9400f/DAR-9500f Storage Media Application Profile Conformance Statement

DAR-9400f/DAR-9500f Media Storage AE conforms to following application profiles.

**Table 19 Supported Application Profile**

| Supported APS       | Real World Activity | Role | SC Option   |
|---------------------|---------------------|------|-------------|
| DAR-9400f/DAR-9500f | Read CD-R           | FSR  | Interchange |
|                     | Write CD-R          | FSC  | Interchange |

#### 3.2.1. Real World Activity - Read CD

The DAR-9400f/DAR-9500f acts as a DICOM FSR with Interchange Service Class Option for images of SOP class in **Table 20**.

**Table 20 Supported SOP classes as FSR**

| SOP Class                        | SOP Class UID                |
|----------------------------------|------------------------------|
| X-ray Angiographic Image Storage | 1.2.840.10008.5.1.4.1.1.12.1 |
| Secondary Capture Image storage  | 1.2.840.10008.5.1.4.1.1.7    |

#### 3.2.1.1. Media Storage Application Profile

Read CD-R applies the following profile.

**Table 21 Supported Application Profile for Read CD-R**

| Supported APS       | Real World Activity | Role | SC Option   |
|---------------------|---------------------|------|-------------|
| DAR-9400f/DAR-9500f | Read CD-R           | FSR  | Interchange |

#### 3.2.1.2. Real World Activity - Write CD

The DAR-9400f/DAR-9500f acts as a DICOM FSC with Interchange Service Class Option for images of SOP class in **Table 22**.

**Table 22 Supported SOP classes as FSR**

| SOP Class                        | SOP Class UID                |
|----------------------------------|------------------------------|
| X-ray Angiographic Image Storage | 1.2.840.10008.5.1.4.1.1.12.1 |
| Secondary Capture Image storage  | 1.2.840.10008.5.1.4.1.1.7    |

#### 3.2.1.3. Media Storage Application Profile

Write CD-R applies the following profile

**Table 23 Supported Application Profile for Write CD-R**

| Supported APS       | Real World Activity | Role | SC Option   |
|---------------------|---------------------|------|-------------|
| DAR-9400f/DAR-9500f | Write CD-R          | FSC  | Interchange |

## 4. Communication Profiles

*DAR-9400f/DAR-9500f* Gateway provides DICOM V3.0 TCP/IP Network Communication Support as defined in Part 8 of the DICOM standard.

### 4.1. TCP/IP Stack

*DAR-9400f/DAR-9500f* Gateway inherits its TCP/IP stack from the computer upon which it executes.

#### 4.1.1. Physical Media Support

*DAR-9400f/DAR-9500f* Gateway is indifferent to the physical medium over which TCP/IP executes; it inherits the medium from the system upon which it executes.

### 4.2. Extensions/Specialization/Privatization

NA

### 4.3. Configuration

#### 4.3.1. AE Title/Presentation Address Mapping

*DAR-9400f/DAR-9500f* maps the Application Entity titles to host name and port number via an internal configuration method. The mapping can be accessed in the configuration menu under the Database tab. Only a privileged user can change the mapping.

*DAR-9400f/DAR-9500f* is connected to network via a router. There are 2 or 3 computers running behind the router, and each computer could be as Storage SCU and Storage Commitment SCU. There is also a computer for Gateway and it could be as Storage SCU.

Therefore, the same IP addresses, up to 4 AE titles and ports settings might be needed for the system which implements the service class SCP.

For other service class, only one computer could be SCU, so an IP address and AE title should be set.

#### 4.3.2. Configurable Parameters

*DAR-9400f/DAR-9500f* Gateway receives its configuration parameters from the user through the AE's GUI.

Configurable parameters are:

1. Local/remote application entity title
2. Local/remote host IP address
3. Local/remote TCP/IP port
4. MAX PDU size
5. Time out for association
6. Route number of study, series, and imageUID
7. Destination of transfer (It is possible for each image, and this is only for Gateway)
8. Switching the output of live image and sub image in DSA (This is only for Gateway)

#### 4.4. Support for Extended Character Sets

*DAR-9400f/DAR-9500f* is known to support the following character sets:

**Table 24 Supported Character Sets**

|                      |                      |
|----------------------|----------------------|
| ISO-IR 100 (default) | Latin Alphabet No. 1 |
| ISO-IR 6             | Basic G0 Set         |
| ISO-IR 87            | Japanese             |

## 5. UID Generation

This section will describe how UID are generated by the *DAR-9400f/DAR-9500f* Gateway system. Gateway has specification of UID generation in different mode, and only one of them can be used at the same time.

### 5.1. Types of Generated Images

The system generates the following images.

#### **X-ray Angiographic Image<sup>1</sup>**

- |                    |   |
|--------------------|---|
| DA image           | : X-ray radiography image (Live image) which is used for cardiac catheterization examination. |
| 3D-DA,3D-DSA image | : X-ray image (Live image) acquired by rotating the C-arm for 3D reconstruction.              |
| DSA image:         | : Subtraction image (Live image) which is used for such as head region study <sup>2</sup> .   |
| Reference image    | : Reference image that represents a selected 1 frame of X-ray radiography image.              |

#### **Secondary Capture Image**

- |            |  |
|------------|--|
| Photo File | : Still image file which is captured from the X-ray radiography image. |
|------------|--|

\* For DA and DSA, those images would be Single-plane or Bi-plane image depending on the system configuration and radiography program.

#### **RDSR**

- |       |  |
|-------|--|
| RDSR: | : Used for study record information such as dose information determined in DICOM standard as Structured Report format. |
|-------|--|

---

<sup>1</sup>Depending on the system configuration and radiography program, images will be either single-plane or bi-plane image for DA and DSA image.

<sup>2</sup>This is a live image normally, however, subtracted image can be sent when using Gateway.

## 5.2. Basic Specification

This section describes the specification of *DAR-9400f/DAR-9500f* without Gateway.

### 5.2.1. Definitions

- Serial Number:** A thirteen digit maximum, number unique to this type of system (*DAR-9400f/DAR-9500f*).
- Study Date:** Date in format YYYYMMDD at which the study was created.
- Study Time:** Time in format HHMMSS at which the study was created.
- Instance Date:** Date in format YYYYMMDD at which the image was created.
- Instance Time:** Time formatted HHMMSS.
- Instance MS:** The milliseconds portion of the time at which the image was created in format mmm.
- Series Number<sup>3</sup>:** Type of the encoding/object:
- 1: DA image (Lossless)
  - 2: DA image (Lossy Cine)
  - 13: Photo File
  - 15: Reference image,
  - 51: DSA image (Lossless Cine)
  - 53: DSA image (Lossy Cine)
- Instance Number<sup>4</sup>:** Sequential Number of the DICOM object generated by the *DAR-9400f/DAR-9500f* for all objects of the same type in the same study.

---

<sup>3</sup>Correspond to 0020,0011 (Series Number).

<sup>4</sup>Correspond to 0020,0013 (Instance Number).

### 5.2.2. Root and Implementation Class UID

*DAR-9400f/DAR-9500f* root is 1.2.392.200036.9110

Acquisition station: Implementation Class UID = <Root>.13.<Serial Number>

Review station: Implementation Class UID = <Root>.14.<Serial Number>

For the Anonymize function:

This root is used only when competitor's files are being anonymized. Otherwise, the above roots are used.

*DAR-9400f/DAR-9500f* Root for anonymization = <Root>.66

Acquisition station: Implementation Class UID = <Root>.66.13.<Serial Number>

Review station: Implementation Class UID = <Root>.66.14.<Serial Number>

### 5.2.3. Study UID<sup>5</sup>

**Study Instance UID** = <ImplementationClassUID>.<StudyDate>.<StudyTime>

### 5.2.4. Series UID

**Series Instance UID** = <StudyInstanceUID>.<SeriesNumber>

### 5.2.5. SOP Instance UID

The Instance sequential number is a number that is generated by the *DAR-9400f/DAR-9500f* sequentially for each new sequence of the same type is a study.

**1- DA, DSA, 3D-DA, 3D-DSA images**

<SeriesInstanceUID>.<Instance Number>

**2- Photo File**

<ImplementationClassUID>.<InstanceDate>.<InstanceTime>.<SeriesNumber>.<InstanceMS>

**3- Reference image**

<ImplementationClassUID>.<InstanceDate>.<InstanceTime>.<SeriesNumber>.<InstanceMS>

**4- Anonymized image**

<ImplementationClassUID(anonymize)>.<SerialNumber>.<InstanceDate>.<InstanceTime>.<SeriesNumber>.<InstanceMS>

---

<sup>5</sup>Use acquired UID when acquiring via MWM.

### 5.3. Renumber Series by Shot Order Mode

This section describes the rules for generating UID when operating in this mode with Gateway.

#### 5.3.1. Series Number

| Image                  | Series Number  |
|------------------------|--|
| DA, DSA, 3D-DSA images | Frontal image: Instance Number *2-1<br>Lateral image: Instance Number *2 |
| Photo File             | 2013   |
| Reference image        | 2015   |

#### 5.3.2. Series Instance UID

- Change the Series Number as in the table above during UID of original image.
- Change the “System Number” to “System Number. 9”.

### 5.4. Unified Series Number Mode

This section describes the rules for generating UID when operating in this mode with Gateway.

#### 5.4.1. Series Number

Change as the following table.

| Image   | Series Number          |
|---------|------------------------|
| DA, DSA | “1” as uniform number. |
| SC      | “1” as uniform number. |
| Ref     | “1” as uniform number. |

#### 5.4.2. Series Instance UID

- Set as same value as the Study Instance UID.
- Change the “System Number” to “System Number. 9”.

## 5.5. Change Series Instance UID Mode

This section describes the rules for generating UID when operating in this mode with Gateway.  
(This mode is related to Photo File transfer. Basic specification is applied to the other images.)

### 5.5.1. Series Instance UID

- Generate the Series Instance UID for each plane and shot.
- To be specific, acquire the value of (0008,1155) Referenced SOP Instance UID in (0008,1140):Referenced Image Sequence of the original image, and overwrite the value “System Number” which is changed to “System Number.9” to the Series Instance UID.

## 5.6. Sort by Shot Order Mode

This section describes the rules for generating UID when operating in this mode with Gateway.  
(This mode is related to Photo File transfer. Basic specification is applied to other images.)

### 5.6.1. Series Number

- Set the Series Number as follows for the Instance Number N of original image.

Frontal:  $2N-1+100$

Lateral:  $2N+100$

### 5.6.2. Instance Number

- Set the Instance Number to (0008,1160) Reference Frame Number.

### 5.6.3. Series Instance UID

- Set the “System Number” to “System Number. 9”
- Change the Series Number as the table above.



## 6. Information Object Definitions

DICOM tags in each image are listed as follows. (This conformance is described on *DAR-9400f Ver.4.2 and DAR-9500f Ver.5.2.40*)

### 6.1. X-ray Angiographic Image

#### 6.1.1. Dicom Meta Information

| Tag         | VR | Type | Attribute Name                     | Remarks |
|-------------|----|------|------------------------------------|---------|
| (0002,0000) | UL | 1    | File Meta Information Group Length |         |
| (0002,0001) | OB | 1    | File Meta Information Version      |         |
| (0002,0002) | UI | 1    | Media Storage SOP Class UID        |         |
| (0002,0003) | UI | 1    | Media Storage SOP Instance UID     |         |
| (0002,0010) | UI | 1    | Transfer Syntax UID                |         |
| (0002,0012) | UI | 1    | Implementation Class UID           |         |
| (0002,0013) | SH | 3    | Implementation Version Name        |         |
| (0002,0016) | AE | 3    | Source Application Entity Title    |         |

#### 6.1.2. Patient Module

| Tag         | VR | Type | Attribute Name       | Remarks |
|-------------|----|------|----------------------|---------|
| (0010,0010) | PN | 2    | Patient's Name       |         |
| (0010,0020) | LO | 2    | Patient ID           |         |
| (0010,0030) | DA | 2    | Patient's Birth Date |         |
| (0010,0040) | CS | 2    | Patient's Sex        |         |

#### 6.1.3. General Study Module

| Tag         | VR | Type | Attribute Name             | Remarks |
|-------------|----|------|----------------------------|---------|
| (0020,000D) | UI | 1    | Study Instance UID         |         |
| (0008,0020) | DA | 2    | Study Date                 |         |
| (0008,0030) | TM | 2    | Study Time                 |         |
| (0008,0050) | SH | 2    | Accession Number           |         |
| (0008,0090) | PN | 2    | Referring Physician's Name |         |
| (0020,0010) | SH | 2    | Study ID                   |         |
| (0008,1030) | LO | 3    | Study Description          |         |

#### 6.1.4. Patient Study Module

| Tag         | VR | Type | Attribute Name   | Remarks |
|-------------|----|------|------------------|---------|
| (0010,1020) | DS | 3    | Patient's Size   |         |
| (0010,1030) | DS | 3    | Patient's Weight |         |

#### 6.1.5. General Series Module

| Tag         | VR | Type | Attribute Name              | Remarks |
|-------------|----|------|-----------------------------|---------|
| (0008,0060) | CS | 1    | Modality                    |         |
| (0020,000E) | UI | 1    | Series Instance UID         |         |
| (0020,0011) | IS | 2    | Series Number               |         |
| (0020,0060) | CS | 2C   | Laterality                  |         |
| (0008,0021) | DA | 3    | Series Date                 |         |
| (0008,0031) | TM | 3    | Series Time                 |         |
| (0008,103E) | LO | 3    | Series Description          |         |
| (0008,1050) | PN | 3    | Performing Physician's Name |         |
| (0008,1070) | PN | 3    | Operator's Name             |         |
| (0018,0015) | CS | 3    | Body Part Examined          |         |
| (0018,5100) | CS | 2C   | Patient Position            |         |

#### 6.1.6. General Equipment Module

| Tag         | VR | Type | Attribute Name                | Remarks |
|-------------|----|------|-------------------------------|---------|
| (0008,0070) | LO | 2    | Manufacturer                  |         |
| (0008,0080) | LO | 3    | Institution Name              |         |
| (0008,0081) | ST | 3    | Institution Address           |         |
| (0008,1010) | SH | 3    | Station Name                  |         |
| (0008,1040) | LO | 3    | Institutional Department Name |         |
| (0018,1020) | LO | 3    | Software Version(s)           |         |

#### 6.1.7. Contrast/Bolus Module

| Tag         | VR | Type | Attribute Name            | Remarks |
|-------------|----|------|---------------------------|---------|
| (0018,0010) | LO | 2    | Contrast/Bolus Agent      |         |
| (0018,1042) | TM | 3    | Contrast/Bolus Start Time |         |

#### 6.1.8. Cine Module

| Tag         | VR | Type | Attribute Name                 | Remarks                           |
|-------------|----|------|--------------------------------|-----------------------------------|
| (0018,1065) | DS | 1C   | Frame Time Vector              | Not available on Reference image. |
| (0008,2144) | IS | 3    | Recommended Display Frame Rate | Not available on Reference image. |
| (0018,0040) | IS | 3    | Cine Rate                      | Only Fluoro image                 |

### 6.1.9. Multi-Frame Module

| Tag         | VR | Type | Attribute Name          | Remarks                           |
|-------------|----|------|-------------------------|-----------------------------------|
| (0028,0008) | IS | 1    | Number of Frames        | Not available on Reference image. |
| (0028,0009) | AT | 1    | Frame Increment Pointer | Not available on Reference image. |

### 6.1.10. Frame Pointers Module

| Tag         | VR | Type | Attribute Name                   | Remarks                           |
|-------------|----|------|----------------------------------|-----------------------------------|
| (0028,6010) | US | 3    | Representative Frame Number      | Not available on Reference image. |
| (0028,6020) | US | 3    | Frame Numbers of Interest (FOI)  | Not available on Reference image. |
| (0028,6022) | LO | 3    | Frame(s) of Interest Description | Not available on Reference image. |

### 6.1.11. Mask Module

(This module is available only if the image is DSA.)

| Tag          | VR | Type | Attribute Name             | Remarks |
|--------------|----|------|----------------------------|---------|
| (0028,6100)  | SQ | 1    | Mask Subtraction Sequence  |         |
| >(0028,6101) | CS | 1    | Mask Operation             |         |
| >(0028,6110) | US | 1C   | Mask Frame Numbers         |         |
| >(0028,6112) | US | 3    | Contrast Frame Averaging   |         |
| >(0028,6114) | FL | 3    | Mask Sub-pixel Shift       |         |
| >(0028,6190) | ST | 3    | Mask Operation Explanation |         |
| (0028,1090)  | CS | 2    | Recommended Viewing Mode   |         |

### 6.1.12. Display Shutter Module

| Tag         | VR | Type | Attribute Name                | Remarks |
|-------------|----|------|-------------------------------|---------|
| (0018,1600) | CS | 1    | Shutter Shape                 |         |
| (0018,1602) | IS | 1C   | Shutter Left Vertical Edge    |         |
| (0018,1604) | IS | 1C   | Shutter Right Vertical Edge   |         |
| (0018,1606) | IS | 1C   | Shutter Upper Horizontal Edge |         |
| (0018,1608) | IS | 1C   | Shutter Lower Horizontal Edge |         |
| (0018,1622) | US | 3    | Shutter Presentation Value    |         |

### 6.1.13. General Image Module

| Tag         | VR | Type | Attribute Name      | Remarks |
|-------------|----|------|---------------------|---------|
| (0020,0013) | IS | 2    | Instance Number     |         |
| (0008,0023) | DA | 2C   | Content Date        |         |
| (0008,0033) | TM | 2C   | Content Time        |         |
| (0020,0020) | CS | 2C   | Patient Orientation |         |
| (0008,0022) | DA | 3    | Acquisition Date    |         |
| (0008,0032) | TM | 3    | Acquisition Time    |         |
| (0020,4000) | LT | 3    | Image Comments      |         |

### 6.1.14. Image Pixel Module

| Tag         | VR | Type | Attribute Name | Remarks |
|-------------|----|------|----------------|---------|
| (0028,0010) | US | 1    | Rows           |         |
| (0028,0011) | US | 1    | Columns        |         |
| (7FE0,0010) | OW | 1    | Pixel Data     |         |

### 6.1.15. X-ray Image Module

| Tag          | VR | Type | Attribute Name               | Remarks   |
|--------------|----|------|------------------------------|---|
| (0008,0008)  | CS | 1    | Image Type                   |   |
| (0028,0002)  | US | 1    | Samples per Pixel            |   |
| (0028,0004)  | CS | 1    | Photometric Interpretation   |   |
| (0028,0100)  | US | 1    | Bits Allocated               |   |
| (0028,0101)  | US | 1    | Bits Stored                  |   |
| (0028,0102)  | US | 1    | High Bit                     |   |
| (0028,0103)  | US | 1    | Pixel Representation         |   |
| (0028,1040)  | CS | 1    | Pixel Intensity Relationship |   |
| (0008,1140)  | SQ | 1C   | Referenced Image Sequence    | Available on Reference, Bi-plane DA and Bi-plane DSA image. |
| >(0008,1150) | UI | 3    | Referenced SOP Class UID     | Available on Reference, Bi-plane DA and Bi-plane DSA image. |
| >(0008,1155) | UI | 3    | Referenced SOP Instance UID  | Available on Reference, Bi-plane DA and Bi-plane DSA image. |

**6.1.16. Curve Module**

| Tag         | VR | Type | Attribute Name            | Remarks                          |
|-------------|----|------|---------------------------|----------------------------------|
| (5000,0005) | US | 1    | Curve Dimensions          | Available when ECG is connected. |
| (5000,0010) | US | 1    | Number of Points          | Available when ECG is connected. |
| (5000,0020) | CS | 1    | Type of Data              | Available when ECG is connected. |
| (5000,0030) | SH | 3    | Axis Units                | Available when ECG is connected. |
| (5000,0103) | US | 1    | Data Value Representation | Available when ECG is connected. |
| (5000,3000) | OW | 1    | Curve Data                | Available when ECG is connected. |
| (5000,0110) | US | 1C   | Curve Data Descriptor     | Available when ECG is connected. |
| (5000,0112) | US | 1C   | Coordinate Start Value    | Available when ECG is connected. |
| (5000,0114) | US | 1C   | Coordinate Step Value     | Available when ECG is connected. |
| (5000,0022) | LO | 3    | Curve Description         | Available when ECG is connected. |
| (5000,0030) | SH | 3    | Axis Units                | Available when ECG is connected. |
| (5000,2500) | LO | 3    | Curve Label               | Available when ECG is connected. |

### 6.1.17. X-Ray Acquisition Module

| Tag         | VR | Type | Attribute Name                        | Remarks  |
|-------------|----|------|---------------------------------------|--|
| (0018,1155) | CS | 1    | Radiation Setting                     |  |
| (0018,0060) | DS | 2    | KVP                                   |  |
| (0018,1147) | CS | 3    | Field of View Shape                   |  |
| (0018,1149) | IS | 3    | Field of View Dimension(s)            |  |
| (0018,1150) | IS | 2C   | Exposure Time                         |  |
| (0018,1151) | IS | 2C   | X-Ray Tube Current                    |  |
| (0018,1152) | IS | 2C   | Exposure                              |  |
| (0018,1154) | DS | 3    | Average Pulse Width                   |  |
| (0018,1155) | CS | 1    | Radiation Setting                     |  |
| (0018,115A) | CS | 3    | Radiation Mode                        |  |
| (0018,115E) | DS | 3    | Image Area Dose Product               | Available when dosimeter is connected, but not available on 3D-DA image. |
| (0018,1164) | DS | 3    | Image Pixel Spacing                   |  |
| (0018,1166) | CS | 3    | Grid                                  |  |
| (0018,1190) | DS | 3    | Focal Spot(s)                         |  |
| (0018,8150) | DS | 3    | Exposure Time (us)                    |  |
| (0028,0030) | DS | 1C   | Pixel Spacing                         |  |
| (0028,0A02) | CS | 3    | Pixel Spacing Calibration Type        |  |
| (0028,0A04) | LO | 1C   | Pixel Spacing Calibration Description |  |
| (0040,8302) | DS | 3    | Entrance Dose in mGy                  |  |

### 6.1.18. X-ray Collimator Module

| Tag         | VR | Type | Attribute Name                       | Remarks |
|-------------|----|------|--------------------------------------|---------|
| (0018,1700) | CS | 1    | Collimator Shape                     |         |
| (0018,1702) | IS | 1C   | Collimator Left Vertical Edge        |         |
| (0018,1704) | IS | 1C   | Collimator Right Vertical Edge       |         |
| (0018,1706) | IS | 1C   | Collimator Upper Horizontal Edge     |         |
| (0018,1708) | IS | 1C   | Collimator Lower Horizontal Edge     |         |
| (0018,1720) | IS | 1C   | Vertices of the Polygonal Collimator |         |

### 6.1.19. XA Positioner Module

| Tag         | VR | Type | Attribute Name                              | Remarks                              |
|-------------|----|------|---|--------------------------------------|
| (0018,1510) | DS | 2    | Positioner Primary Angle                    |                                      |
| (0018,1511) | DS | 2    | Positioner Secondary Angle                  |                                      |
| (0018,1500) | CS | 2C   | Positioner Motion                           |                                      |
| (0018,1520) | DS | 2C   | Positioner Primary Angle Increment          | Available on 3D-DA image.            |
| (0018,1521) | DS | 2C   | Positioner Secondary Angle Increment        | Available on 3D-DA image.            |
| (0018,1110) | DS | 3    | Distance Source to Detector                 |                                      |
| (0018,1111) | DS | 3    | Distance Source to Patient                  |                                      |
| (0018,1114) | DS | 3    | Estimated Radiographic Magnification Factor |                                      |
| (0018,1530) | DS | 3    | Detector Primary Angle                      | Available on 3D and rotational image |
| (0018,1531) | DS | 3    | Detector Secondary Angle                    | Available on 3D and rotational image |

### 6.1.20. SOP Common Module

| Tag         | VR | Type | Attribute Name         | Remarks |
|-------------|----|------|------------------------|---------|
| (0008,0005) | CS | 1C   | Specific Character Set |         |
| (0008,0016) | UI | 1C   | SOP Class UID          |         |
| (0008,0018) | UI | 1C   | SOP Instance UID       |         |
| (0008,0012) | DA | 3    | Instance Creation Date |         |
| (0008,0013) | TM | 3    | Instance Creation Time |         |

### 6.1.21. VOI LUT Module

| Tag         | VR | Type | Attribute Name | Remarks |
|-------------|----|------|----------------|---------|
| (0028,1051) | DS | 1C   | Window Width   |         |
| (0028,1050) | DS | 3    | Window Center  |         |

### 6.1.22. Additional Attributes Module

| Tag         | VR | Type | Attribute Name           | Remarks                       |
|-------------|----|------|--------------------------|-------------------------------|
| (0008,1160) | IS | 3    | Referenced Frame Number  | Available on Reference image. |
| (0018,7050) | CS | 3    | Filter Material          |                               |
| (0018,7052) | DS | 3    | Filter Thickness Minimum |                               |
| (0018,7054) | DS | 3    | Filter Thickness Maximum |                               |

## 6.2. Secondary Capture Image

### 6.2.1. Dicom Meta Information

| Tag         | VR | Type | Attribute Name                     | Remarks |
|-------------|----|------|------------------------------------|---------|
| (0002,0000) | UL | 1    | File Meta Information Group Length |         |
| (0002,0001) | OB | 1    | File Meta Information Version      |         |
| (0002,0002) | UI | 1    | Media Storage SOP Class UID        |         |
| (0002,0003) | UI | 1    | Media Storage SOP Instance UID     |         |
| (0002,0010) | UI | 1    | Transfer Syntax UID                |         |
| (0002,0012) | UI | 1    | Implementation Class UID           |         |
| (0002,0013) | SH | 3    | Implementation Version Name        |         |
| (0002,0016) | AE | 3    | Source Application Entry Title     |         |

### 6.2.2. Patient Module

| Tag         | VR | Type | Attribute Name       | Remarks |
|-------------|----|------|----------------------|---------|
| (0010,0010) | PN | 2    | Patient's Name       |         |
| (0010,0020) | LO | 2    | Patient ID           |         |
| (0010,0030) | DA | 2    | Patient's Birth Date |         |
| (0010,0040) | CS | 2    | Patient's Sex        |         |

### 6.2.3. General Study Module

| Tag         | VR | Type | Attribute Name             | Remarks |
|-------------|----|------|----------------------------|---------|
| (0020,000D) | UI | 1    | Study Instance UID         |         |
| (0008,0020) | DA | 2    | Study Date                 |         |
| (0008,0030) | TM | 2    | Study Time                 |         |
| (0008,0050) | SH | 2    | Accession Number           |         |
| (0008,0090) | PN | 2    | Referring Physician's Name |         |
| (0020,0010) | SH | 2    | Study ID                   |         |
| (0008,1030) | LO | 3    | Study Description          |         |

### 6.2.4. Patient Study Module

| Tag         | VR | Type | Attribute Name   | Remarks |
|-------------|----|------|------------------|---------|
| (0010,1020) | DS | 3    | Patient's Size   |         |
| (0010,1030) | DS | 3    | Patient's Weight |         |



### 6.2.5. General Series Module

| Tag         | VR | Type | Attribute Name              | Remarks |
|-------------|----|------|-----------------------------|---------|
| (0020,000E) | UI | 1    | Series Instance UID         |         |
| (0020,0011) | IS | 2    | Series Number               |         |
| (0020,0060) | CS | 2C   | Laterality                  |         |
| (0008,0021) | DA | 3    | Series Date                 |         |
| (0008,0031) | TM | 3    | Series Time                 |         |
| (0008,103E) | LO | 3    | Series Description          |         |
| (0008,1050) | PN | 3    | Performing Physician's Name |         |
| (0008,1070) | PN | 3    | Operators' Name             |         |
| (0018,0015) | CS | 3    | Body Part Examined          |         |
| (0018,5100) | CS | 2C   | Patient Position            |         |

### 6.2.6. General Equipment Module

| Tag         | VR | Type | Attribute Name                | Remarks |
|-------------|----|------|-------------------------------|---------|
| (0008,0070) | LO | 2    | Manufacturer                  |         |
| (0008,0080) | LO | 3    | Institution Name              |         |
| (0008,0081) | ST | 3    | Institution Address           |         |
| (0008,1010) | SH | 3    | Station Name                  |         |
| (0008,1040) | LO | 3    | Institutional Department Name |         |
| (0018,1020) | LO | 3    | Software Version(s)           |         |

### 6.2.7. SC Equipment Module

| Tag         | VR | Type | Attribute Name  | Remarks |
|-------------|----|------|-----------------|---------|
| (0008,0064) | CS | 1    | Conversion Type |         |
| (0008,0060) | CS | 3    | Modality        |         |

### 6.2.8. General Image Module

| Tag          | VR | Type | Attribute Name              | Remarks |
|--------------|----|------|-----------------------------|---------|
| (0020,0013)  | IS | 2    | Instance Number             |         |
| (0008,0023)  | DA | 2C   | Content Date                |         |
| (0008,0033)  | TM | 2C   | Content Time                |         |
| (0020,0020)  | CS | 2C   | Patient Orientation         |         |
| (0008,0008)  | CS | 3    | Image Type                  |         |
| (0008,0022)  | DA | 3    | Acquisition Date            |         |
| (0008,0032)  | TM | 3    | Acquisition Time            |         |
| (0008,2111)  | ST | 3    | Derivation Description      |         |
| (0020,4000)  | LT | 3    | Image Comments              |         |
| (0008,1140)  | SQ | 3    | Referenced Image Sequence   |         |
| >(0008,1150) | UI | 3    | Referenced SOP Class UID    |         |
| >(0008,1155) | UI | 3    | Referenced SOP Instance UID |         |

### 6.2.9. Image Pixel Module

| Tag         | VR | Type | Attribute Name             | Remarks |
|-------------|----|------|----------------------------|---------|
| (0028,0002) | US | 1    | Samples per Pixel          |         |
| (0028,0004) | CS | 1    | Photometric Interpretation |         |
| (0028,0010) | US | 1    | Rows                       |         |
| (0028,0011) | US | 1    | Columns                    |         |
| (0028,0100) | US | 1    | Bits Allocated             |         |
| (0028,0101) | US | 1    | Bits Stored                |         |
| (0028,0102) | US | 1    | High Bit                   |         |
| (0028,0103) | US | 1    | Pixel Representation       |         |
| (7FE0,0010) | OW | 1    | Pixel Data                 |         |

### 6.2.10. SC Image Module

| Tag         | VR | Type | Attribute Name                        | Remarks |
|-------------|----|------|---------------------------------------|---------|
| (0028,0030) | DS | 1C   | Pixel Spacing                         |         |
| (0028,0A02) | CS | 3    | Pixel Spacing Calibration Type        |         |
| (0028,0A04) | LO | 1C   | Pixel Spacing Calibration Description |         |

### 6.2.11. X-ray Collimator Module

| Tag         | VR | Type | Attribute Name                       | Remarks |
|-------------|----|------|--------------------------------------|---------|
| (0018,1700) | CS | 1    | Collimator Shape                     |         |
| (0018,1720) | IS | 1C   | Vertices of the Polygonal Collimator |         |

### 6.2.12. SOP Common Module

| Tag         | VR | Type | Attribute Name         | Remarks |
|-------------|----|------|------------------------|---------|
| (0008,0005) | CS | 1C   | Specific Character Set |         |
| (0008,0016) | UI | 1C   | SOP Class UID          |         |
| (0008,0018) | UI | 1C   | SOP Instance UID       |         |
| (0008,0012) | DA | 3    | Instance Creation Date |         |
| (0008,0013) | TM | 3    | Instance Creation Time |         |

### 6.2.13. VOI LUT Module

| Tag         | VR | Type | Attribute Name | Remarks |
|-------------|----|------|----------------|---------|
| (0028,1051) | DS | 1C   | Window Width   |         |
| (0028,1050) | DS | 3    | Window Center  |         |

#### 6.2.14. Additional Attributes Module

| Tag         | VR | Type | Attribute Name               | Remarks |
|-------------|----|------|------------------------------|---------|
| (0008,1160) | IS | 3    | Referenced Frame Number      |         |
| (0018,1154) | DS | 3    | Average Pulse Width          |         |
| (0018,1160) | SH | 3    | Filter Type                  |         |
| (0018,7050) | CS | 3    | Filter Material              |         |
| (0018,7052) | DS | 3    | Filter Thickness Minimum     |         |
| (0018,7054) | DS | 3    | Filter Thickness Maximum     |         |
| (0028,1040) | CS | 3    | Pixel Intensity Relationship |         |
| (0040,8302) | DS | 3    | Entrance Dose in mGy         |         |

### 6.3. RDSR

#### 6.3.1. Dicom Meta Information

| Tag         | VR | Type | Attribute Name                     | Remarks |
|-------------|----|------|------------------------------------|---------|
| (0002,0000) | UL | 1    | File Meta Information Group Length |         |
| (0002,0001) | OB | 1    | File Meta Information Version      |         |
| (0002,0002) | UI | 1    | Media Storage SOP Class UID        |         |
| (0002,0003) | UI | 1    | Media Storage SOP Instance UID     |         |
| (0002,0010) | UI | 1    | Transfer Syntax UID                |         |
| (0002,0012) | UI | 1    | Implementation Class UID           |         |
| (0002,0013) | SH | 3    | Implementation Version Name        |         |
| (0002,0016) | AE | 3    | Source Application Entity Title    |         |

#### 6.3.2. Patient Module

| Tag         | VR | Type | Attribute Name       | Remarks |
|-------------|----|------|----------------------|---------|
| (0010,0010) | PN | 2    | Patient's Name       |         |
| (0010,0020) | LO | 2    | Patient ID           |         |
| (0010,0030) | DA | 2    | Patient's Birth Date |         |
| (0010,0040) | CS | 2    | Patient's Sex        |         |

#### 6.3.3. General Study Module

| Tag         | VR | Type | Attribute Name             | Remarks |
|-------------|----|------|----------------------------|---------|
| (0020,000D) | UI | 1    | Study Instance UID         |         |
| (0008,0020) | DA | 2    | Study Date                 |         |
| (0008,0030) | TM | 2    | Study Time                 |         |
| (0008,0050) | SH | 2    | Accession Number           |         |
| (0008,0090) | PN | 2    | Referring Physician's Name |         |
| (0020,0010) | SH | 2    | Study ID                   |         |
| (0008,1030) | LO | 3    | Study Description          |         |

#### 6.3.4. Patient Study Module

| Tag         | VR | Type | Attribute Name   | Remarks |
|-------------|----|------|------------------|---------|
| (0010,1020) | DS | 3    | Patient's Size   |         |
| (0010,1030) | DS | 3    | Patient's Weight |         |

### 6.3.5. SR Document Series Module

| Tag         | VR | Type | Attribute Name      | Remarks |
|-------------|----|------|---------------------|---------|
| (0008,0060) | CS | 1    | Modality            |         |
| (0020,000E) | UI | 1    | Series Instance UID |         |
| (0020,0011) | IS | 2    | Series Number       |         |
| (0008,0021) | DA | 3    | Series Date         |         |
| (0008,0031) | TM | 3    | Series Time         |         |
| (0008,103E) | LO | 3    | Series Description  |         |

### 6.3.6. General Equipment Module

| Tag         | VR | Type | Attribute Name                | Remarks |
|-------------|----|------|-------------------------------|---------|
| (0008,0080) | LO | 3    | Institution Name              |         |
| (0008,0081) | ST | 3    | Institution Address           |         |
| (0008,1010) | SH | 3    | Station Name                  |         |
| (0008,1040) | LO | 3    | Institutional Department Name |         |

### 6.3.7. Enhanced General Equipment Module

| Tag         | VR | Type | Attribute Name            | Remarks |
|-------------|----|------|---------------------------|---------|
| (0008,0070) | LO | 2    | Manufacturer              |         |
| (0008,1090) | LO | 3    | Manufacturer's Model Name |         |
| (0018,1000) | LO | 3    | Device Serial Number      |         |
| (0018,1020) | LO | 3    | Software Version(s)       |         |

### 6.3.8. SR Document General Module

| Tag         | VR | Type | Attribute Name    | Remarks |
|-------------|----|------|-------------------|---------|
| (0020,0013) | IS | 2    | Instance Number   |         |
| (0040,A491) | CS | 1    | Completion Flag   |         |
| (0040,A493) | CS | 1    | Verification Flag |         |
| (0008,0023) | DA | 2C   | Content Date      |         |
| (0008,0033) | TM | 2C   | Content Time      |         |

### 6.3.9. SOP Common Module

| Tag         | VR | Type | Attribute Name         | Remarks |
|-------------|----|------|------------------------|---------|
| (0008,0005) | CS | 1C   | Specific Character Set |         |
| (0008,0016) | UI | 1    | SOP Class UID          |         |
| (0008,0018) | UI | 1    | SOP Instance UID       |         |

### 6.3.10. Additional Attributes Module

| Tag         | VR | Type | Attribute Name              | Remarks |
|-------------|----|------|-----------------------------|---------|
| (0008,1050) | PN | 3    | Performing Physician's Name |         |
| (0008,1070) | PN | 3    | Operators' Name             |         |
| (0018,5100) | CS | 3    | Patient Position            |         |
| (0040,A050) | CS | 3    | Continuity of Content       |         |
| (0040,A504) | SQ | 3    | Content Template Sequence   |         |
| (0008,0105) | CS | 3    | Mapping Resource            |         |
| (0040,DB00) | CS | 3    | Template Identifier         |         |

### 6.3.11. SR Document Content Module

#### 6.3.11.1. TID 10001 Projection X-Ray Radiation Dose

| VT        | Concept Name                                   | VM  | Value Set Constraint                                    | Remarks                                    |
|-----------|--|-----|---|--|
| CONTAINER | EV(113701, DCM, "X-Ray Radiation Dose Report") | 1   |   |  |
| CODE      | EV(121058, DCM, "Procedure reported")          | 1   | DT(113704, DCM, "Projection X-Ray")                     |  |
| CODE      | EV(G-C0E8, SRT, "Has Intent")                  | 1   | DCID(3629) procedure Intent                             |  |
| INCLUDE   | DTID(1002) Observer Context                    | 1-n |   | Refer to DTID(1002).                       |
| CODE      | EV(113705, DCM, "Scope of Accumulation")       | 1   | DCID(10000) Scope of Accumulation                       |  |
| UIDREF    | DCID(10001) UID Types                          | 1   |   |  |
| INCLUDE   | DTID(10002) Accumulated X-Ray Dose             | 1   | EV(113622, DCM, "Single Plane")                         | Refer to DTID(10002). For Single-plane     |
| INCLUDE   | DTID(10002) Accumulated X-Ray Dose             | 1   | EV(113620, DCM, "Plane A")                              | Refer to DTID(10002). For Bi-plane Frontal |
| INCLUDE   | DTID(10002) Accumulated X-Ray Dose             | 1   | EV(113621, DCM, "Plane B")                              | Refer to DTID(10002). For Bi-plane Lateral |
| INCLUDE   | DTID(10003) Irradiation Event X-Ray Data       | 1-n |   |  |
| CODE      | EV(113854, DCM, "Source of Dose Information")  | 1-n | DCID(10020) Source of Projection X-Ray Dose Information |  |

**6.3.11.2. TID 1002 Observer Context**

| VT      | Concept Name   | VM | Value Set Constraint    | Remarks              |
|---------|--|----|-------------------------|----------------------|
| CODE    | EV(121005, DCM, "Observer Type")                     | 1  | DCID(270) Observer Type |                      |
| INCLUDE | DTID(1004)<br>Device observer identifying attributes | 1  |                         | Refer to DTID (1004) |

**6.3.11.3. TID 1004 Device Observer Identifying Attributes**

| VT     | Concept Name                                     | VM | Value Set Constraint   | Remarks                  |
|--------|--|----|--|--------------------------|
| UIDREF | EV(121012, DCM, "Device Observer UID")           | 1  |  | Implementation Class UID |
| TEXT   | EV(121013, DCM, "Device Observer Name")          | 1  | Defaults to value of Station Name(0008,1010) in General Equipment Module             |                          |
| TEXT   | EV(121014, DCM, "Device Observer Manufacturer")  | 1  | Defaults to value of Manufacturer (0008,0070) in General Equipment Module            |                          |
| TEXT   | EV(121015, DCM, "Device Observer Model Name")    | 1  | Defaults to value of Manufacturer's Model Name(0008,1090) in GeneralEquipment Module |                          |
| TEXT   | EV (121016,DCM, "Device Observer Serial Number") | 1  | Defaults to value of DeviceSerial Number (0018,1000) in General Equipment Module     |                          |



**6.3.11.4. TID 10002 Accumulated X-Ray Dose**

| VT        | Concept Name                                     | VM  | Value Set Constraint                 | Remarks   |
|-----------|--|-----|--------------------------------------|---|
| CONTAINER | EV(113702, DCM, "Accumulated X-Ray Dose Data")   | 1   |                                      |   |
| INCLUDE   | EV(113764, DCM, "Acquisition Plane")             | 1   |                                      | Either one of the following:<br>113622, DCM, Single Plane<br>113620, DCM, "PlaneA"<br>113621, DCM, "PlaneB" |
| CONTAINER | EV(122505, DCM, "Calibration")                   | 1-n |                                      |   |
| CODE      | EV(113794,DCM, "Dose measurement")               | 1   | DCID(10010) Dose measurement Devices |   |
| DATETIME  | EV(113723, DCM, "Calibration Date")              | 1   |                                      |   |
| NUM       | EV(122322, DCM, "Calibration Factor")            | 1   | Units = EV(1, UCUM, "no units")      |   |
| NUM       | EV(113763, DCM, "Calibration Uncertainty")       | 1   | Units = EV(%, UCUM, "Percent")       |   |
| TEXT      | EV(113724, DCM, "Calibration Responsible Party") | 1   |                                      |   |
| INCLUDE   | DTID(10004) Accumulated Projection X-Ray Dose    | 1   |                                      | Refer to DTID(10004).   |

**6.3.11.5. TID 10004 Accumulated Projection X-Ray Dose**

| VT  | Concept Name   | VM  | Value Set Constraint             | Remarks  |
|-----|--|-----|----------------------------------|--|
| NUM | EV(113722, DCM, "Dose Area Product Total")             | 1   | Units = EV(Gym2, UCUM, "Gym2")   |  |
| NUM | EV(113725, DCM, "Dose (RP) Total")                     | 1   | Units = EV(Gy, UCUM, "Gy")       |  |
| NUM | EV(113726, DCM, "Fluoro Dose Area Product Total")      | 1   | Units = EV(Gy.m2, UCUM, "Gy.m2") |  |
| NUM | EV(113728, DCM, "Fluoro Dose (RP) Total")              | 1   | Units = EV(Gy, UCUM, "Gy")       |  |
| NUM | EV(113730, DCM, "Total Fluoro Time")                   | 1   | Units = EV(s, UCUM, "s")         |  |
| NUM | EV(113727, DCM, "Acquisition Dose Area product Total") | 1   | Units = EV(Gy.m2, UCUM, "Gy.m2") |  |
| NUM | EV(113729, DCM, "Acquisition Dose (RP) Total")         | 1   | Units = EV(Gy, UCUM, "Gy")       |  |
| NUM | EV(113855, DCM, "Total Acquisition Time")              | 1   | Units = EV(s, UCUM, "s")         |  |
| NUM | EV(113731, DCM, "Total Number of Radiographic Frames") | 1   | Units = EV(1, UCUM, "no units")  |  |
| NUM | DCID(10008)<br>Dose related Distance Measurements      | 1-n | Units = EV(mm, UCUM, "mm")       | Set [Distance Source to Reference Point] of the first or the last acquisition. |

**6.3.11.6. TID 10003 Irradiation Event X-Ray Data**

| VT        | Concept Name                                    | VM | Value Set Constraint                          | Remarks   |
|-----------|---|----|---|---|
| CONTAINER | EV(113706, DCM, "Irradiation Event X-Ray Data") | 1  |   |   |
| CODE      | EV(113764, DCM, "Acquisition Plane ")           | 1  | DCID(10003)<br>Equipment Plane Identification | Either one of the following:<br>113620, DCM, Plane A<br>113621, DCM, Plane B<br>113622, DCM, Single |
| DATETIME  | DT(111526, DCM, "DataTime Started")             | 1  |   |   |
| CODE      | EV(113721, DCM, "Irradiation Event Type")       | 1  | DCID(10002)<br>Irradiation Event Type         |   |
| TEXT      | EV(125203, DCM, "Acquisition Protocol")         | 1  |   | Set DUP name during acquisition.  |
| CODE      | EV(113780, DCM, "Reference Point Definition")   | 1  |   |   |
| UIDREF    | EV(113769, DCM, "Irradiation Event UID")        | 1  |   |   |
| NUM       | EV(122130, DCM, "Dose Area Product")            | 1  | Units = EV(Gy.m2, UCUM, "Gy.m2")              |   |
| NUM       | EV(113738, DCM, "Dose(RP)")                     | 1  | Units = EV(Gy, UCUM, "Gy")                    |   |
| NUM       | EV(112011, DCM, Positioner Primary Angle)       | 1  | Units = EV(deg, UCUM, "deg")                  |   |
| NUM       | EV(112012, DCM, "Positioner Secondary Angle")   | 1  | Units = EV(deg, UCUM, deg)                    |   |

| VT        | Concept Name                                      | VM  | Value Set Constraint                   | Remarks                         |
|-----------|---|-----|--|---------------------------------|
| NUM       | EV(113739, DCM, "Positioner Primary End Angle")   | 1   | Units = EV(deg, UCUM, "deg")           |                                 |
| NUM       | EV(113740, DCM, "Positioner Secondary End Angle") | 1   | Units = EV(deg, UCUM, "deg")           |                                 |
| NUM       | EV(113790, DCM, "Collimated Field Area")          | 1   | Units = EV(m2, UCUM, "m^2")            |                                 |
| CONTAINER | EV(113771, DCM, "X-Ray Filter")                   | 1-n |  |                                 |
| CODE      | EV(113772, DCM, "X-Ray Filter Type")              | 1   | DCID(10007) X-Ray Filter Types         |                                 |
| CODE      | EV(113757, DCM, "X-Ray Filter material")          | 1   | DCID(10006) X-Ray Filter Material      |                                 |
| NUM       | EV(113758, DCM, "X-Ray Filter Thickness Minimum") | 1   | Units = EV(mm, UCUM, "mm")             |                                 |
| NUM       | EV(113773, DCM, "X-Ray Filter Thickness Maximum") | 1   | Units = EV(mm, UCUM, "mm")             |                                 |
| CODE      | EV(113732, DCM, "Fluoro Mode")                    | 1   | DCID(10004) Fluoro Modes               |                                 |
| NUM       | EV(113791, DCM, "Pulse Rate")                     | 1   | Units = EV({pulse}/s, UCUM, "pulse/s") |                                 |
| NUM       | EV(113768, DCM, "Number of Pulses")               | 1   | Units = EV(1, UCUM, "no units")        |                                 |
| NUM       | EV(113733, DCM "KVP")                             | 1-n | EV(kV, UCUM, "kV")                     | Configure only for radiography. |
| NUM       | EV(113724, DCM, "X-Ray Tube Current")             | 1-n | Units = EV(ms, UCUM, "ms")             | Configure only for radiography. |
| NUM       | EV(113824, DCM, "Exposure Time")                  | 1-n | Units = EV(ms, UCUM, "ms")             | Configure Only for radiography. |
| NUM       | EV(113793, DCM, "Pulse Width")                    | 1-n | Units = EV(ms, UCUM, "ms")             |                                 |
| NUM       | EV(113736, DCM, Exposure)                         | 1-n | Units = EV(uAs, UCUM, "uAs")           |                                 |
| NUM       | EV(113766, DCM, "Focal Spot Size")                | 1   | Units = EV(mm, UCUM, "mm")             |                                 |
| NUM       | EV(113742, DCM, "Irradiation Duration")           | 1   | Units = EV(s, UCUM, "s")               |                                 |
| CODE      | EV(113745, DCM, "Patient table Relationship")     | 1   | DCID(21) Patient Gantry Relationship   |                                 |

| VT    | Concept Name                                      | VM  | Value Set Constraint                  | Remarks |
|-------|---|-----|---------------------------------------|---------|
| CODE  | EV(113743, DCM, "Patient Orientation")            | 1   | DCID(19) Patient Orientation          |         |
| CODE  | EV(113744, DCM, "Patient Orientation Modifier")   | 1   | DCID(20) Patient Orientation Modifier |         |
| NUM   | DCID(10008)<br>Dose Related Distance Measurements | 1-n | Units = EV(mm, UCUM, "mm")            |         |
| CODE  | EV(123014, DCM, "Target Region")                  | 1   | DCID(4031) Common Anatomic Region     |         |
| IMAGE | EV(113795, DCM, "Acquired Image")                 | 1-n |                                       |         |

## 6.4. Data Dictionary of Private Attributes

All private attributes are deleted from images by executing the image processing on the gateway system.

| Tag          | VR | Type | Attribute Name                 | Remarks |
|--------------|----|------|--------------------------------|---------|
| (0029,0015)  | LO | 1    | Private Creator                |         |
| (0029,1002)  | DS | —    | Private Time Vector            |         |
| (0029, 1501) | DS | 1    | Tilting Angle                  |         |
| (0029, 1502) | IS | 1    | FPD Size                       |         |
| (0029, 1509) | LO | 1    | DUP Name                       |         |
| (0029,1517)  | US | 1    | Pre-Processing Horizontal Flip |         |
| (0029,1518)  | US | 1    | Pre-Processing Vertical Flip   |         |
| (7FDF,0010)  | LO | —    | Header Padding Group           |         |
| (7FDF,1001)  | OB | —    | Header Padding                 |         |

Use the following tags only for radiography images.

| Tag         | VR | Type | Attribute Name | Remarks |
|-------------|----|------|----------------|---------|
| (0029,1516) | ID | —    | Rad Type       |         |

Use the following tags for StentView images.

| Tag         | VR | Type | Attribute Name           | Remarks |
|-------------|----|------|--------------------------|---------|
| (0029,1513) | LO | 1    | StentView                |         |
| (0029,1530) | US | 1    | Number of StentView ROIs |         |
| (0029,1531) | US | 1C   | Rows of Proximal ROI     |         |
| (0029,1532) | US | 1C   | Number of StentView ROIs |         |
| (0029,1533) | LO | 1C   | Position of Proximal ROI |         |
| (0029,1534) | US | 1C   | Rows of Distal ROI       |         |
| (0029,1535) | US | 1C   | Columns of Distal ROI    |         |
| (0029,1536) | LO | 1C   | Position of Distal ROI   |         |

Use the following tags for 3D images.

| Tag          | VR | Type | Attribute Name  | Remarks |
|--------------|----|------|---|---------|
| (0029,1506)  | DS | 1    | Field Of View   |         |
| (0029,1507)  | DS | 1    | Distance of Source to Detector                        |         |
| (0029,1508)  | DS | 1    | Distance of Source to Patient                         |         |
| (0029,1528)  | DS | 1    | Table Top Vertical Position                           |         |
| (0029,1529)  | DS | 1    | Table Top Longitudinal Position                       |         |
| (0029,152A)  | DS | 1    | Table Top Lateral Position                            |         |
| (0029,152D)  | DS | 1    | Ceiling Travel Longitudinal Position                  |         |
| (0029,152E)  | DS | 1    | Ceiling Travel Transversal Position                   |         |
| (0029,152F)  | DS | 1    | ISO Center Height                                     |         |
| (6B01, 0001) | LO | —    | Private Creator                                       |         |
| (6B01, 0100) | LO | —    | 3D-DSA, 3D-DA, 3D-RSM, (CB Future)                    |         |
| (6B01, 0101) | CS | —    | Sensor Type: 0=I.I, 1=FPD                             |         |
| (6B01, 0102) | CS | —    | 08:MH200, 09:MH300, 10:MH200S                         |         |
| (6B01, 0103) | DS | —    | Mask Frames Count/Frames Before Rotation/Total Frames |         |
| (6B01, 0108) | DS | —    | Rotation Speed (degrees/sec)                          |         |
| (6B01,010A)  | IS | 3    | BH Filter   |         |
| (6B01,0110)  | IS | —    | 3D Reconstruction Mode                                |         |
| (6B01, 0180) | OB | —    | mAs   |         |